

Lims1 Web Portal

Improve Customer Retention Free Up Your Staff

Lims1 Web Portal provides your customers instant access to their Laboratory Data

- Laboratory Data is available 24/7. Customers can access at any time, from anywhere
- Clients can retrieve historical results, Certificates of Analysis. The Laboratory Database effectively becomes your customer's database
- With the Lims1 Sample Submission Add-On, Clients can Register Samples prior to Delivery. This can speed up delivery, reduce lab staff data entry, reduce errors.

Completely Rebuilt from Scratch


Lims1 Web Portal replaces KIWI Web Portal. It has been completely rebuilt from the ground up using the latest Technology and Tools. Your customers will see a difference and your IT Department will love the improvements.

The screenshot shows the Lims1 Web Portal interface. On the left is a dark navigation sidebar with the Lims1 logo and user information (Vivek User). The main area features a table with columns for Lab Ref, Description, Due Date, Reg Date, Final Date, Client Ref, Status, Contact, and Download Pdf. The table contains several rows of data related to National Audit Testing.

Lab Ref	Description	Due Date	Reg Date	Final Date	Client Ref	Status	Contact	Download Pdf
T1500066	National Audit Paper Kraft	23/06/2016	09/06/2016		BLANK	Incomplete		[Download]
T1500063	Board National Audit Testing	23/06/2016	09/06/2016		BLANK	Complete		[Download]
T1500062	Board National Audit Testing	23/06/2016	09/06/2016		BLANK	Incomplete		[Download]
T1500061	Board National Audit Testing	23/06/2016	09/06/2016		BLANK	Complete		[Download]
T1500060	Board National Audit Testing	23/06/2016	09/06/2016		BLANK	Incomplete		[Download]
T1500058	Board National Audit Testing	22/06/2016	08/06/2016		BLANK	No Entries		[Download]
T1500057	Board National Audit Testing	22/06/2016	08/06/2016		BLANK	Incomplete		[Download]

Intuitive User Friendly Screens

- Modern Presentation
- Highly Customisable
- Very Secure (SSL)
- Latest Technology
- Mobile Ready
- Registration of Samples



New and Improved Features

Many improvements have been incorporated into the new version in response to requests from our existing customers.

Easy Assignment of User Roles and Privileges

- Can easily administer Users, Roles and Privileges from Browser based Admin.
- Clients can add and Edit users
- Users can manage their own passwords

Logout Button

Logout button now accessible from all screens for added security

Mobile Friendly Screens

Use paragraphs often. A wall of white text makes it hard for the reader to skim a story and find a way to quickly drop in and out of your content. White space gives the user an opening into your information. Don't be afraid to leave spaces open.

Remote Sample Registration

An optional Feature is now available for remote registration of samples and testing requirements. Customers can pre register samples to greatly enhance communication, registration and resource planning.

NOTE: Requires Lims1 Sample Submission Add-On

Mobile Ready

Lims1 Web Portal is ready for mobile phone access. The portal automatically adjusts for the user device.

Clients can now easily access their data in the field.

With the inclusion of the optional Sample Registration feature clients can log in samples as soon as they have collected them in the field. This reduces paperwork for the client and the lab. The lab has prior notice of expected samples and the sample registration is significantly streamlined.



Mobile Ready



Lims1 - Built with Proven Technologies

The Technical Bits

Development Platform	Server Side: Microsoft Latest Technologies Client Side: AngularJS, Bootstrap, JQuery, HTML, Responsive for Mobile Compatibility
Deployment Platform	IIS 7.5 or above
Hardware Requirements	Server which can host IIS 7.5 or above
Security	Role and permission based authentication and authorization SSL available
Compliance	Built to comply with OWASP guidelines

Recent Implementations

International Exposure

AgriGen's scientific expertise in PCR method development and validation and their 'regulatory awareness' deliver state-of-the-art, world class DNA testing services.

AgriGen customers are located in all corners of the globe. The Web Portal provides the ability for their clients to access results and files at any time of the day or night from any time zone.

Contact Us

Give us a call for more information about our services and products

LTech Australia
5/7 Packard Avenue
Castle Hill NSW 2154

+612 8677 9565

sales@lms1.com

Visit us on the web at
www.lms1.com

Lims1 Web Portal

Providing laboratories with easy & flexible
LIMS products and services, for superior lab management.

Detailed Features and Availability

August, 2016

Feature	Comment	Availability
Users and Login	Username and password	Now
	Can Limit Menus by user	Now
	Can restrict data by client Reference	Now
	Log out button	Now
	Add Users	Now
	Modify Users	Now
	Suspend/Delete Users	Now
	User can modify password	Now
	Configurable Password Policy	Now
	Assign Users to Roles	Now
Branding	Can modify logos, styles, images	Now – included in Installation
Can Associate with different Data-bases	Required for DEV, UAT and PROD	Now – multiple instances
Column Properties	Change Order	Now – requires LTech assistance
	Sort	Now
	Search	Now
	Change Headings	Now – requires LTech assistance
	Change Column Source	Now – requires LTech assistance
Links to Documents	Links to pdf if exists	Now
	Links to other docs if exists	Now – included in Installation
Appearance	Modern Interface	Now
	Mobile Ready	Now
	Highlighted Cells	Now – requires LTech assistance
	Good Response time	Now
Graphing	Ability to Graph Data	Coming Soon
Export	Export to Excel	Coming Soon
Register Samples	Register or Preregister	Now – Requires Lims1 API